DOSEUR

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DESCRIPTION

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INSTRUCTIONS FOR USE

DESCRIPTION

Measuring cups are used for measurement and administration of syrups. Measuring cups are made by UNION PLASTIC by injection of plastic materials consistent for food contact and/ or European pharmacopoeia. Some references may be colored by the use of a coloring master batch food contact compliant. The cups are used according to the dosage indicated by the pharmaceutical manufacturer on the package leaflet attached the drug. Read this leaflet before use. Measuring cups can be single or multiple uses for a single patient. The device provided should be used with the particular drug with which it is included.

USE

Before the first use, it is necessary to wash manually the cup with water.

The cups must be used in an environment allowing the reading of the graduations.

1. **Wash hands before use**
2. **Dosage:**
 - Detach the measuring cup from the stopper of the medicine bottle,
 - Check that the measuring cup is not damaged and well cleaned,
 - Open the bottle,
 - Put syrups drug into the associated cups to the mark corresponding to the indicated dosage.
3. **Administration:**
 - Release the contents of the cup into the patient's mouth.
4. **Cleaning after use:**
 - Wash manually the cup with water,
 - Rinse with water until the remaining medical residues disappear,
 - Dry the cup,
 - Store the cup with its syrup in its box in a place inaccessible to children.

INDICATIONS FOR USE

The cups should be used according to the procedure of the pharmaceutical manufacturer and according to the drug to be administered. This medical device is reusable until the end of treatment. The instructions for use should be translated into the official languages of the country of use in the package leaflet. Any serious incident involving the device should be notified to the manufacturer and to the competent authority of the Member State in which the user and / or patient is established.

CONTRAINDICATIONS

- Do not use the cup if it has been badly cleaned or damaged.
- Do not use the cup for another application or for any medicine other than the one used.
- Do not use for administration in other ports.
- The cup cannot be used beyond the expiry date of the drug.

DEGREE OF MEDICAL DEVICE ACCURACY

The general accuracy of the cups is $\pm 10\%$ for all volumes.

PRECAUTIONS OF USE

DO NOT USE THE DEVICE IF IT IS DAMAGED.

The washing temperature must not exceed 50°C: do not wash in the dishwasher, do not boil. The product should not be used by people who have no insanity of mind or impaired mental faculties and can read and understand the instructions. Keep out of children. Check that the packaging is not damaged.

CONDITIONS OF STORAGE

The cup should be kept clean and dry in the medicine's packaging at a temperature between 5°C and 50°C. The cup must be disposed of at the same time as the container of the drug with household waste, or return to the pharmacist all the opened boxes including the cup as well as the vials for an appropriate destruction.

MANUFACTURER INFORMATION

First placing on the market: 1998
Date of last revision of the instructions for use: Notice DT2 revision D 2020/10/29
Number of the notified body for the CE marking identification: 0459

UNION PLASTIC SAS
ZA Robert, 43140 St-Didier en Velay – France | Phone: +33 (0)4 71 61 13 09 | union-plastic@omerin.com

MEASURING CUP

INSTRUCTIONS FOR USE

DESCRIPTION

Measuring cups are used for measurement and administration of syrups.

Measuring cups are made by UNION PLASTIC by injection of plastic materials consistent for food contact and/ or European pharmacopeia.

Some references may be colored by the use of a coloring master batch food contact compliant.

The cups are used according to the dosage indicated by the pharmaceutical manufacturer on the package leaflet attached the drug.

Read this leaflet before use. Measuring cups can be single or multiple uses for a single patient.

The device provided should be used with the particular drug with which it is included.

USE

Before the first use, it is necessary to wash manually the cup with water.

The cups must be used in an environment allowing the reading of the graduations.

1. Wash hands before use

2. Dosage:

- Detach the measuring cup from the stopper of the medicine bottle,
- Check that the measuring cup is not damaged and well cleaned,
- Open the bottle,
- Put syrups drug into the associated cups to the mark corresponding to the indicated dosage.

3. Administration:

- Release the contents of the cup into the patient's mouth.

4. Cleaning after use:

- Wash manually the cup with water,
- Rinse with water until the remaining medical residues disappear,
- Dry the cup,
- Store the cup with its syrup in its box in a place inaccessible to children.

INDICATIONS FOR USE

The cups should be used according to the procedure of the pharmaceutical manufacturer and according to the drug to be administered.

This medical device is reusable until the end of treatment.

The instructions for use should be translated into the official languages of the country of use in the package leaflet.

Any serious incident involving the device should be notified to the manufacturer and to the competent authority of the Member State in which the user and / or patient is established.

CONTRAINDICATIONS

- Do not use the cup if it has been badly cleaned or damaged.
- Do not use the cup for another application or for any medicine other than the one used.
- Do not use for administration in other ports.
- The cup cannot be used beyond the expiry date of the drug.

DEGREE OF MEDICAL DEVICE ACCURACY

The general accuracy of the cups is $\pm 10\%$ for all volumes.

PRECAUTIONS OF USE

DO NOT USE THE DEVICE IF IT IS DAMAGED.

The washing temperature must not exceed 50°C: do not wash in the dishwasher, do not boil.

The product should not be used by people who have no insanity of mind or impaired mental faculties and can read and understand the instructions.

Keep out of children.

Check that the packaging is not damaged.

CONDITIONS OF STORAGE

The cup should be kept clean and dry in the medicine's packaging at a temperature between 5°C and 50°C.

The cup must be disposed of at the same time as the container of the drug with household waste, or return to the pharmacist all the opened boxes including the cup as well as the vials for an appropriate destruction.

MANUFACTURER INFORMATION

First placing on the market: 1998

Date of last revision of the instructions for use: Notice DT2 revision D 2020/10/29

Number of the notified body for the CE marking identification: 0459

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INSTRUCTIONS FOR USE

DESCRIPTION

Measuring cups are used for measurement and administration of syrups.
Measuring cups are made by UNION PLASTIC by injection of plastic materials consistent for food contact and/ or European pharmacopeia.
Some references may be colored by the use of a coloring master batch food contact compliant.
The cups are used according to the dosage indicated by the pharmaceutical manufacturer on the package leaflet attached the drug.
Read this leaflet before use. Measuring cups can be single or multiple uses for a single patient.
The device provided should be used with the particular drug with which it is included.

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INDICATIONS FOR USE

The cups should be used according to the procedure of the pharmaceutical manufacturer and according to the drug to be administered.
This medical device is reusable until the end of treatment.
The instructions for use should be translated into the official languages of the country of use in the package leaflet.
Any serious incident involving the device should be notified to the manufacturer and to the competent authority of the Member State in which the user and / or patient is established.

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The washing temperature must not exceed 50°C: do not wash in the dishwasher, do not boil.
The product should not be used by people who have no insanity of mind or impaired mental faculties and can read and understand the instructions.
Keep out of children.
Check that the packaging is not damaged.

CONDITIONS OF STORAGE

The cup should be kept clean and dry in the medicine's packaging at a temperature between 5°C and 50°C.
The cup must be disposed of at the same time as the container of the drug with household waste, or return to the pharmacist all the opened boxes including the cup as well as the vials for an appropriate destruction.

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USE

Before the first use, it is necessary to wash manually the cup with water.
The cups must be used in an environment allowing the reading of the graduations.

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PRECAUTIONS OF USE

DO NOT USE THE DEVICE IF IT IS DAMAGED.
The washing temperature must not exceed 50°C: do not wash in the dishwasher, do not boil.
The product should not be used by people who have no insanity of mind or impaired mental faculties and can read and understand the instructions.
Keep out of children.
Check that the packaging is not damaged.

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1. **Wash hands before use**
2. **Dosage:**
 - Detach the measuring cup from the stopper of the medicine bottle,
 - Check that the measuring cup is not damaged and well cleaned,
 - Open the bottle,
 - Put syrups drug into the associated cups to the mark corresponding to the indicated dosage.
3. **Administration:**
 - Release the contents of the cup into the patient' s mouth.
4. **Cleaning after use:**
 - Wash manually the cup with water,
 - Rinse with water until the remaining medical residues disappear,
 - Dry the cup,
 - Store the cup with its syrup in its box in a place inaccessible to children.

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